

# Diffraction of Light Through a Plane Grating, Single Slit, Wire and Circular Aperture Lab Report 2

Ayaan Dutt  
Lab partner: Raheem Manoj

Professor Bikram Phookun  
TAs: Ankit Shrestha, Chandra Shekhar Mishra, Kamal Nayan  
TF: Uzma Khan

Date of Experiment: February 20th, 2023

Date of Submission: February 27th, 2023

## Aim

- Part A: To determine the wavelength  $\lambda$  of the light emitted by a laser source by studying the diffraction of light due to plane diffraction gratings.
- Part B: To determine the width of the given single slit by studying its diffraction pattern.
- Part C: To determine the diameter of a given wire/hair strand by studying its diffraction pattern.
- Part D: To determine the size of the circular aperture by studying its diffraction pattern.
- Part E: To determine the thickness, pitch and angle of a single helical spring by studying its diffraction pattern.

## Experimental Setup

List of apparatus:

- Optical bench
- 10 mW semiconductor red laser source with a mount
- 5 mW DPSS green laser source with a mount
- Measuring tape
- Small plastic scale with least count 0.5 mm
- Graph sheet
- Masking tape
- Set of plane diffraction gratings with grating spacing 100 lines/mm, 200 lines/mm and 600 lines/mm with a holder and mounts
- Single slit of fixed width mounted on a slide
- Thin wire/hair strand mounted on a slide

- Set of circular apertures mounted on a slide
- Single helix (spring) set in a holder with a mount
- Spirit level

Least count of measuring tape = 0.001 m

Least count of small plastic scale = 0.0005 m




Figure 1: Experimental setup for the diffraction of light

## Theoretical Background

### Single Slit (or Thin Wire)

According to Babinet's Principle, complementary objects produce the same diffraction pattern, except for the intensity of the central maxima. This implies that a narrow single slit and a single thin wire/hair strand are complementary since they produce the same diffraction pattern and can be understood by the same theory.

Single slit diffraction (also known as, Fraunhofer diffraction) gives a pattern of varying intensity consisting of a bright central maximum with alternate minima and maxima of decreasing intensity on either side. This is caused by the constructive and destructive interference of light waves passing through the slit.




Figure 2: Schematic diagram of diffraction through a single slit (Fraunhofer diffraction)

The position of minimas is given by the equation:

$$d \sin \theta_m = \pm m \lambda \text{ for } m = 1, 2, 3 \dots \quad (1)$$

Where  $d$  is the width of the slit,  $\theta_m$  is the angle corresponding to the  $m$ th minima and  $\lambda$  is the wavelength of the incident light. The  $\pm$  indicates which side of the central maxima the minima lies on. For small angles we can approximate  $\sin \theta_m \approx \tan \theta_m = x_m/D$  where  $x_m$  is the distance between the central maxima and the  $m$ th minima on the screen and  $D$  is the distance between the slit and the screen.

## Plane Diffraction Grating

A transmission diffraction grating consists of a large number of slits separated from one another by an opaque region. The diffraction pattern of a plane diffraction grating consists of widely separated distinct bright intense spots known as principal maxima with faint secondary maxima scattered between them.




Figure 3: Schematic diagram of diffraction through a plane transmission grating

The location of principal maxima on the screen is given by:

$$d \sin \theta_m = m\lambda \text{ for } m = 1, 2, 3 \dots \quad (2)$$

Where  $d$  is the slit width (calculated as  $1/\text{number of lines per meter}$ ),  $\theta_m$  is the angle corresponding to the  $m$ th principal maxima and  $\lambda$  is the wavelength of the incident light. Here too, we can approximate  $\sin \theta_m \approx \tan \theta_m = x_m/D$  where  $x_m$  is the distance between the central maxima and the  $m$ th minima on the screen and  $D$  is the distance between the slit and the screen. Hence, we can modify our equation as:

$$\frac{dx_m}{D} = m\lambda \implies x_m = \frac{m\lambda D}{d} \text{ for } m = 1, 2, 3 \dots \quad (3)$$

## Circular Aperture




Figure 4: Schematic diagram of diffraction through a circular aperture

The diffraction pattern produced by a circular aperture is known as Airy diffraction, which is similar to Fraunhofer diffraction but is mathematically more complex. However, the same expression can be used with a slight change in the coefficient:

$$d \sin \theta_m = \bar{m} \lambda \text{ for } \bar{m} = 1.22, 2.23, 3.23, 4.24 \dots \quad (4)$$

Where  $d$  is the diameter of the aperture,  $\theta_m$  is the angle corresponding to the  $m$ th dark ring and  $\lambda$  is the wavelength of the incident light.

## Single Helical Spring




Figure 5: (a) Typical view of helical spring  
(b) Schematic diagram when viewed at normal incidence<sup>1</sup>

A helical spring has 3 measurable parameters - thickness of wire ( $a_1$ ), distance between loops ( $d_1$ ) and angle  $2\alpha_1$ . When viewed at normal incidence its projection is equivalent to two sets of parallel wires of the same thickness at an angle  $2\alpha_1$  to each other.




Figure 6: Diffraction pattern produced by a single helical spring<sup>2</sup>

Each line of the X-shaped diffraction pattern corresponds to one set of parallel wires where the broader set of minima are due diffraction from a single wire while the narrow set of minima are formed by the interference due to multiple wires. The angle between the two arms is equal to  $2\alpha_1$ .

We can use equation (3) to find the distance between two consecutive minima:

$$\beta = x_{m+1} - x_m = \frac{(m+1)\lambda D}{d} - \frac{m\lambda D}{d}$$

$$\Rightarrow \beta = \frac{\lambda D}{d}$$

The two length scales in the diffraction pattern can be associated with two fringe widths -  $\beta_1$  and  $\beta_2$ , such that:

$$a_1 = \frac{\lambda D}{\beta_1} \text{ and } d_1 = \frac{\lambda D}{\beta_2} \quad (5)$$

The pitch of a spring can be found from the distance between loops  $d_1$  and the angle of the spring  $2\alpha_1$  using:

$$p = \frac{d_1}{\cos \alpha_1} \quad (6)$$

<sup>1</sup>Praveen Pathak, Charudatt Kadolkar, and Manish Kapoor. *Diffraction due to Helical Structure*. 2015. URL: <https://ipho.olimpicos.net/pdf/IPhO2015Q4.pdf>

<sup>2</sup>Physics Alive. "Episode 43 - DNA Physics and DNA of Class." Physics Alive!, 31 May 2022, URL: [physic-salive.com/dna](https://physic-salive.com/dna).

## Procedure

**Step 1:** Set up the laser mount and diffraction slide mount on the optical bench such that they are at the same height and aligned horizontally. Check this with the spirit level.

**Step 2:** Tape a graph sheet on the screen at the central maxima to accurately mark and measure the diffraction pattern. Measure the distance between the diffraction slide and the screen with a tape measure.

*[Precaution: The distance between the diffraction slide and screen should be significantly more than the width of the diffraction object (in the order of  $10^4$ - $10^6$  times larger) to be able to distinguish a clear diffraction pattern.]*

**Step 3:** Using a plane diffraction grating with known slit width, measure the distance between principal maxima in the diffraction pattern to determine the wavelength  $\lambda$  of the light source using equation (2).

*[Precaution: For proper functioning, keep the laser on throughout the experiment. When not in use, use a light-blocking screen to block the beam.]*

**Step 4:** Using a single slit diffraction slide, measure the distance between minima in the diffraction pattern produced. Calculate the slit width using equation (1).

**Step 5:** Using a mounted thin wire/hair strand, measure the distance between minima in the diffraction pattern produced. Calculate the thickness of the wire/hair using equation (1).

**Step 6:** Using a circular aperture diffraction slide, measure the diameter of the dark rings in the diffraction pattern produced. Calculate the aperture diameter using equation (4).

**Step 7:** Using a mounted single helical spring, measure the primary and secondary fringe widths in the diffraction pattern produced. Also measure the angle between the arms of the X-shaped pattern. Calculate the thickness of the wire, pitch and angle of the spring using equation (5).

## Part A: Plane Diffraction Grating

### Observations

Grating Number ( $\text{mm}^{-1}$ )	Screen Distance $D$ (m)	Slit Width $d$ (m)	1st Maxima Distance $x_1$ (m)	2nd Maxima Distance $x_2$ (m)	3rd Maxima Distance $x_3$ (m)
100	2.283	$1 \times 10^{-5}$	0.124	0.252	0.375
200	2.335	$5 \times 10^{-6}$	0.248	0.514	0.794
600	2.328	$1.667 \times 10^{-6}$	0.775	1.972	-

Table 1: The distance to principal maxima in the diffraction pattern of green light through a plane grating

To find the approximate wavelength of red light, we took only one reading:

Grating Number =  $200 \text{ mm}^{-1}$

Screen Distance ( $D$ ) = 2.324 m

Slit Width ( $d$ ) =  $5 \times 10^{-6} \text{ m}$

1st Maxima Distance ( $x_1$ ) = 0.315 m

## Analysis

We can rearrange the values in equation (2) to find the wavelength of light:

$$\lambda = \frac{d \sin \theta}{m} \text{ for } m = 1, 2, 3 \dots$$

Where,

$$\theta = \arctan\left(\frac{x_m}{D}\right)$$

Grating No. (mm <sup>-1</sup> )	Theta (rad)			Wavelength $\lambda$ (m)			Average Wavelength $\lambda_{avg}$ (m)
	$\theta_1$	$\theta_2$	$\theta_3$	$\lambda_1$	$\lambda_2$	$\lambda_3$	
100	0.0543	0.1078	0.1614	$5.427 \cdot 10^{-7}$	$5.381 \cdot 10^{-7}$	$5.356 \cdot 10^{-7}$	$5.388 \cdot 10^{-7}$
200	0.1082	0.2167	0.3277	$5.399 \cdot 10^{-7}$	$5.374 \cdot 10^{-7}$	$5.365 \cdot 10^{-7}$	$5.379 \cdot 10^{-7}$
600	0.3272	0.7027	-	$5.358 \cdot 10^{-7}$	$5.387 \cdot 10^{-7}$	-	$5.373 \cdot 10^{-7}$

Table 2: Wavelength of green light calculated from the diffraction pattern of a plane grating

Average wavelength of green light ( $\lambda_g$ ) =  $5.380 \cdot 10^{-7}$  m

For red light,

$$\theta = 0.1347$$

$$\lambda_r = \frac{d \sin \theta}{m} = 6.715 \cdot 10^{-7} \text{ m}$$

## Part B: Single Slit

### Observations

Screen Distance ( $D$ ) = 2.374 m

Wavelength $\lambda$ (m)	Distance between 1st minima $y_1$ (cm)	Distance between 2nd minima $y_2$ (cm)
$5.380 \cdot 10^{-7}$	2.7	5.6
$6.715 \cdot 10^{-7}$	3.6	7.0

Table 3: Distance between minima for single slit diffraction with red and green laser light

### Analysis

The width of the slit was found using the equation (1):

$$d = \frac{\pm m \lambda}{\sin \theta_m} \text{ for } m = 1, 2, 3 \dots$$

Where  $\theta_m = \tan^{-1}(x_m/D)$  and  $x_m = y_m/2$  since  $y_m$  is the distance between corresponding minima.



Wavelength $\lambda$ (m)	Slit Width $d_1$ (mm)	Slit Width $d_2$ (mm)	Avg Slit Width $d_{avg}$ (mm)
$5.380 \cdot 10^{-7}$	0.0870	0.0908	0.0889
$6.715 \cdot 10^{-7}$	0.0888	0.0913	0.09005

Table 4: Slit width calculated from Fraunhofer diffraction pattern

Average Slit Width ( $d$ ) = 0.0894 mm

## Part C: Thin Wire/Hair Strand

### Observations

Screen Distance ( $D$ ) = 2.374 m

Wavelength $\lambda$ (m)	Distance between 1st minima $y_1$ (cm)	Distance between 2nd minima $y_2$ (cm)
$5.380 \cdot 10^{-7}$	2.9	7.2
$6.715 \cdot 10^{-7}$	3.2	7.0

Table 5: Distance between minima for diffraction of a hair strand with red and green laser light

### Analysis

Due to Babinet's Principle, we can use the same equation as single slit diffraction for a hair strand. The diameter of the hair was found using equation (1):

$$d = \frac{\pm m \lambda}{\sin \theta_m} \text{ for } m = 1, 2, 3 \dots$$

Where  $\theta_m = \tan^{-1}(x_m/D)$  and  $x_m = y_m/2$  since  $y_m$  is the distance between corresponding minima.

Wavelength $\lambda$ (m)	Hair Diameter $d_1$ (mm)	Hair Diameter $d_2$ (mm)	Avg Hair Diameter $d_{avg}$ (mm)
$5.380 \cdot 10^{-7}$	0.0879	0.0708	0.0794
$6.715 \cdot 10^{-7}$	0.0999	0.0914	0.0957

Table 6: Hair diameter calculated from Fraunhofer diffraction pattern

Average Hair Diameter = 0.0876 mm.

## Part D: Circular Aperture

### Observations

Screen Distance( $D$ ) = 2.374 m

Circle No.	Wavelength $\lambda$ (m)	Diameter of 1st dark ring $y_1$ (cm)	Diameter of 2nd dark ring $y_2$ (cm)
1	$5.380 \cdot 10^{-7}$	0.9	1.6
1	$6.715 \cdot 10^{-7}$	1.1	2.0
2	$5.380 \cdot 10^{-7}$	0.5	1.8
2	$6.715 \cdot 10^{-7}$	0.6	1.0

Table 7: Diameter of dark rings for diffraction through a circular aperture with red and green laser light

### Analysis

The diameter of each circular aperture was determined using equation (4):

$$d = \frac{\bar{m}\lambda}{\sin \theta_m} \text{ for } \bar{m} = 1.22, 2.23, 3.23, 4.24 \dots$$

Where  $\theta_m = \tan^{-1}(x_m/D)$  and  $x_m = y_m/2$  since  $y_m$  is the diameter and  $x_m$  is the radius.

Circle No.	Wavelength $\lambda$ (m)	Diameter $d_1$ (mm)	Diameter $d_2$ (mm)	Avg Diameter $d_{avg}$ (mm)
1	$5.380 \cdot 10^{-7}$	0.355	0.357	0.356
1	$6.715 \cdot 10^{-7}$	0.346	0.356	0.351
2	$5.380 \cdot 10^{-7}$	0.650	0.713	0.682
2	$6.715 \cdot 10^{-7}$	0.622	0.711	0.667

Table 8: Width of aperture 1 using equation (4)

Average diameter of the circular aperture 1 = 0.354 mm

Average diameter of the circular aperture 2 = 0.675 mm

## Part E: Single Helical Spring

### Observations

Wavelength of green light ( $\lambda_g$ ) =  $5.380 \cdot 10^{-7}$

Screen Distance ( $D$ ) = 2.374 m

Primary Fringe Width ( $\beta_1$ ) = 0.85 mm =  $8.5 \cdot 10^{-4}$

Secondary Fringe Width ( $\beta_2$ ) = 0.15 mm =  $1.5 \cdot 10^{-4}$

Angle between diffraction arms ( $2\alpha_1$ ) =  $24^\circ = 0.4189$  rad

## Analysis

Using equations (5), we can determine the thickness of the wire ( $a_1$ ) and the distance between parallel loops ( $d_1$ ):

$$a_1 = \frac{\lambda D}{\beta_1} = 1.503 \cdot 10^{-4} m = 0.1503 mm$$

$$d_1 = \frac{\lambda D}{\beta_2} = 8.515 \cdot 10^{-4} m = 0.8515 mm$$

The pitch of the spring can be determined from equation (6):

$$p = \frac{d_1}{\cos \alpha_1} = 8.690 \cdot 10^{-4} m = 0.869 mm$$

## Error Analysis

We can do error analysis for the value of wavelength  $\lambda$  for a given value of  $d$  if we use equation (3):

$$\lambda = \frac{x_m d}{m D} \text{ for } m = 1, 2, 3 \dots$$

Where we can use the following formula to find the error:

$$\frac{\Delta \lambda}{\lambda} = \frac{\Delta x_m}{x_m} + \frac{\Delta d}{d} + \frac{\Delta D}{D}$$

As there is no error in the measurement of  $d$  since it is calculated from the given grating number, we can eliminate that term:

$$\begin{aligned} \frac{\Delta \lambda}{\lambda} &= \frac{\Delta x_m}{x_m} + \frac{\Delta D}{D} \\ &= \frac{10^{-3}}{0.124} + \frac{10^{-3}}{2.283} \\ &= 8.06 \cdot 10^{-3} + 4.38 \cdot 10^{-4} \\ &= 8.498 \cdot 10^{-3} \end{aligned}$$

$$\Delta \lambda = 4.61 \cdot 10^{-9} m = 4.61 nm$$

Note that the relative error due to  $x_m$  is one order of magnitude higher than the relative error due to  $D$ , which suggests that we must reduce our  $\Delta x_m$  (i.e. use a more precise measuring instrument) to improve our experiment.

## Results

### Part A

$$\text{Wavelength of Green Light } (\lambda_g) = \boxed{5.380 \cdot 10^{-7} m}$$

$$\text{Wavelength of Red Light } (\lambda_r) = \boxed{6.715 \cdot 10^{-7} m}$$

### Part B

$$\text{Width of Single Slit } (a) = \boxed{0.0894 mm}$$

### Part C

Diameter of Hair Strand ( $d$ ) =  $0.0876mm$

### Part D

Diameter of the Circular Aperture 1 ( $d_1$ ) =  $0.354mm$

Diameter of the Circular Aperture 2 ( $d_2$ ) =  $0.675mm$

### Part E

Thickness of Single Helical Spring Wire ( $a_1$ ) =  $0.1503mm$

Pitch of Single Helical Spring ( $p$ ) =  $0.869mm$

Angle of Single Helical Spring ( $\alpha_1$ ) =  $12^\circ = 0.2094rad$

## Discussion

To increase precision in the experiment, the following improvements can be made:

- The entire experimental setup should be put in a dark room so that the diffraction pattern is bright and distinct.
- A digital measuring instrument with a least count of 0.01 mm would enhance the results of the measurements.
- A larger screen to allow for the collection of more data points so that we can graph the data and arrive at a better value for  $\lambda$  and  $d$ .